
py-lmd

Release 1.0.0

Georg Wallmann, Sophia Mädler and Niklas Schmacke

Jan 31, 2022

CONTENTS:

1	Quick Start	2
1.1	Installation from Github	2
1.2	Generating Shapes	2
1.3	Using the py-lmd tools	5
1.4	Numbers and Letters	7
1.5	Text	10
1.6	Different Coordinate Systems	11
2	Modules	13
2.1	lmd.lib	13
	2.1.1 Collection	13
	2.1.2 Shape	14
2.2	lmd.tools	15
	Python Module Index	18
	Index	19

Table of Contents

- *Overview*

Welcome to the show

QUICK START

1.1 Installation from Github

To install the py-lmd library clone the Github repository and use pip to install the library in your current environment. It is recommended to use the library with a conda environment. Please make sure that the package is installed editable like described. Otherwise static glyph files might not be available.

```
git clone https://github.com/HornungLab/py-lmd
pip install -e .
```

Once installed the modules can be loaded as following:

```
from lmd.lib import Collection, Shape
```

1.2 Generating Shapes

As first example we will create a cutting data for a rectangle and visualize it. First we need to think of a calibration points for our coordinate system. The calibration points are specified as numpy array and should have the shape (3, 2). When calibrating the file on the Leica LMD after loading the file, the order of selecting the calibration points is the same as the order of points set here. If no other orientation transform has been specified, the calibration points are defined in the (x, y) coordinate system. More on the usage of different coordinate systems can be found under XXX.

```
import numpy as np
from lmd.lib import Collection, Shape

calibration = np.array([[0, 0],
                        [0, 100],
                        [50, 50]])
```

With these calibration coordinates we can create our *Collection* object. The *Collection* is the base object for creating cutting data and holds a list of *Shape* objects. Furthermore, it allows to read and write the Leica LMD XML format and handles the coordinate system.


```
my_first_collection = Collection(calibration_points = calibration)
```

We can then create our first rectangle by using the *Shape* object and passing it to our collection by using the *add_shape* method. For creating the Shape object we need to pass the vertices as numpy array. List of vertices should always be closed with the last vertex equaling the first one.

```
rectangle_coordinates = np.array([[10,10],
                                  [40,10],
                                  [40,40],
                                  [10,40],
                                  [10,10]])
rectangle = Shape(rectangle_coordinates)
my_first_collection.add_shape(rectangle)
```

We can visualize our collection with the `plot` method. Using the `calibration = True` parameter will include the calibration coordinates in the plot.


```
my_first_collection.plot(calibration = True)
```



We can generate a second shape in form of a triangle. This time we will be using the `new_shape` method of the collection object.

```
triangle_coordinates = np.array([[10,70], [40,70], [40,100], [10,70]])
my_first_collection.new_shape(triangle_coordinates)

my_first_collection.plot(calibration = True)
```



We can then export and save our collection of shapes into xml cutting data.

```
my_first_collection.save("first_collection.xml")
```

```
<?xml version='1.0' encoding='UTF-8'?>
<ImageData>
  <GlobalCoordinates>1</GlobalCoordinates>
  <X_CalibrationPoint_1>0</X_CalibrationPoint_1>
  <Y_CalibrationPoint_1>0</Y_CalibrationPoint_1>
  <X_CalibrationPoint_2>0</X_CalibrationPoint_2>
  <Y_CalibrationPoint_2>10000</Y_CalibrationPoint_2>
  <X_CalibrationPoint_3>5000</X_CalibrationPoint_3>
  <Y_CalibrationPoint_3>5000</Y_CalibrationPoint_3>
  <ShapeCount>2</ShapeCount>
  <Shape_1>
    <PointCount>5</PointCount>
    <X_1>1000</X_1>
    <Y_1>1000</Y_1>
    <X_2>4000</X_2>
    <Y_2>1000</Y_2>
    <X_3>4000</X_3>
    <Y_3>4000</Y_3>
```

(continues on next page)

(continued from previous page)

```

    <X_4>1000</X_4>
    <Y_4>4000</Y_4>
    <X_5>1000</X_5>
    <Y_5>1000</Y_5>
  </Shape_1>
  <Shape_2>
    <PointCount>4</PointCount>
    <X_1>1000</X_1>
    <Y_1>7000</Y_1>
    <X_2>4000</X_2>
    <Y_2>7000</Y_2>
    <X_3>4000</X_3>
    <Y_3>10000</Y_3>
    <X_4>1000</X_4>
    <Y_4>7000</Y_4>
  </Shape_2>
</ImageData>

```

Looking at the generated xml output we can see the calibration points and different shapes. Furthermore, we see that the coordinate system has been scaled by a linear scaling factor. As all points are defined as integers scaling by a linear factor allows to use decimal numbers as coordinates.

1.3 Using the py-lmd tools

A lot of useful functionality is included in the tools module of the py-lmd package. We will first use the rectangle functionality to create rectangle shapes fast.

```

import numpy as np
from lmd.lib import Collection, Shape
from lmd import tools

calibration = np.array([[0, 0], [0, 100], [50, 50]])
my_first_collection = Collection(calibration_points = calibration)


```

After initiating the coordinate system we can use the `rectangle()` helper function to create a `Shape` object with a rectangle with specified size and position.

```


my_square = tools.rectangle(10, 10, offset=(10,10))
my_first_collection.add_shape(my_square)
my_first_collection.plot(calibration = True)

```



We can further specify an angle of rotation.

```
my_square = tools.rectangle(20, 20, offset=(30,30), rotation = np.pi/4)
my_first_collection.add_shape(my_square)
my_first_collection.plot(calibration = True)
```

1.4 Numbers and Letters


The py-lmd tools offer a limited support for numbers and some capital letters. The following glyphs are available: *ABCDEFGHI0123456789-.* They were included in the package as they allow for the development of more consistent calibration and sample indexing. In screens with multiple slides, samples can be unambiguously identified from imaged data.

We will first use `glyphs()` to load single glyphs. The glyphs are included in the py-lmd package as SVG files and are loaded by the `svg_to_lmd()` into an uncalibrated *Collection*. This uncalibrated collection is returned and can be joined with a calibrated collection with the `join()` function.

```
import numpy as np
from lmd.lib import Collection, Shape
from lmd import tools


calibration = np.array([[0, 0], [0, 100], [50, 50]])
my_first_collection = Collection(calibration_points = calibration)

digit_1 = tools.glyph(1)
my_first_collection.join(digit_1)
my_first_collection.plot(calibration = True)
```




By default glyphs and text have a height of ten units and are located by the top left corner. We can use the *offset* and *multiplier* parameters to change the size and position.

```
digit_2 = tools.glyph(2, offset = (0,80), multiplier = 5)
my_first_collection.join(digit_2)
my_first_collection.plot(calibration = True)
```



Like with the previous rectangle example we can also use the *rotation* parameter to set a clockwise rotation.

```
glyph_A = tools.glyph('A', offset=(0,80), rotation =-np.pi/4)
my_first_collection.join(glyph_A)
my_first_collection.plot(calibration = True)
```




1.5 Text

Next to individual glyphs the `text()` method can be used to write text with specified position, size and rotation.

```
import numpy as np
from lmd.lib import Collection, Shape
from lmd import tools


calibration = np.array([[0, 0], [0, 100], [100, 50]])
my_first_collection = Collection(calibration_points = calibration)

identifier_1 = tools.text('0456_B2', offset=np.array([30, 40]), rotation = -np.pi/4)
my_first_collection.join(identifier_1)
my_first_collection.plot(calibration = True)
```



1.6 Different Coordinate Systems

The coordinates for the Leica LMD are defined as (x, y) coordinates with the x-axis extending to the right and the y-axis extending to the top. All cutting data should exist in this coordinate system and should be calibrated accordingly. When cutting data is generated based on whole slide images we have to keep in mind that images are often indexed differently. Images in Fiji or Numpy are indexed as $(row, column)$ with the rows extending downwards and the columns extending to the right. If we want to identify positions in image data - like calibration crosses or single cells - we have to translate their position in the $(row, column)$ format to the (x, y) format.



Conversion

$s_f = 10$ # scale parameter
 $\mathbf{0} = \begin{bmatrix} 0 & -1 \\ 1 & 0 \end{bmatrix}$ # orientation_transform parameter
 $(3, 1) \times \mathbf{0} \times s_f = (10, -30)$

The py-lmd library has been designed in a way which allows to transform the coordinate system prior to saving. Therefore one can specify all coordinates in the image coordinate system and rely on the library to handle the transformation. In this case the *orientation_transform* attribute needs to be set when the Collection is created.

```
calibration = np.array([[10, 10], [1000, 10], [500, 500]])

collection = Collection(calibration_points = calibration)
collection.orientation_transform = np.array([[0, -1], [1, 0]])
```

In this case all coordinates for calibration points and shapes can be set in form of (*row*, *column*) coordinates. The orientation transform is only applied when the Collection is saved or, if desired, when the Collection is plotted.

MODULES

2.1 lmd.lib

2.1.1 Collection

class `lmd.lib.Collection`(*calibration_points: Optional[numpy.ndarray] = None*)

Class which is used for creating shape collections for the Leica LMD6 & 7. Contains a coordinate system defined by calibration points and a collection of various shapes.

Parameters `calibration_points` – Calibration coordinates in the form of (3, 2).

shapes

Contains all shapes which are part of the collection.

Type `List[Shape]`

calibration_points

Calibration coordinates in the form of (3, 2).

Type `Optional[np.ndarray]`

orientation_transform

defines transformations performed on the provided coordinate system prior to export as XML. This orientation_transform is always applied to shapes when there is no individual orientation_transform provided.

Type `np.ndarray`

add_shape(*shape: lmd.lib.Shape*)

Add a new shape to the collection.

Parameters `shape` – Shape which should be added.

join(*collection: lmd.lib.Collection*)

Join the collection with the shapes of a different collection. The calibration markers of the current collection are kept. Please keep in mind that coordinate systems and calibration points must be compatible for correct joining of collections.

Parameters `collection` – Collection which should be joined with the current collection object.

load(*file_location: str*)

Can be used to load a shape file from XML. Both, XMLs generated with py-lmd and the Leica software can be used. :param file_location: File path pointing to the XML file.

new_shape(*points: numpy.ndarray, well: Optional[str] = None, name: Optional[str] = None*)

Directly create a new Shape in the current collection.

Parameters

- **points** – Array or list of lists in the shape of $(N,2)$. Contains the points of the polygon forming a shape.
- **well** – Well in which to sort the shape after cutting. For example A1, A2 or B3.
- **name** – Name of the shape.

plot(*calibration: bool = True, mode: str = 'line', fig_size: tuple = (5, 5), apply_orientation_transform: bool = True, apply_scale: bool = False, save_name: Optional[str] = None*)

This function can be used to plot all shapes of the corresponding shape collection.

Parameters

- **calibration** – Controls whether the calibration points should be plotted as crosshairs. Deactivating the crosshairs will result in the size of the canvas adapting to the shapes. Can be especially useful for small shapes or debugging.
- **fig_size** – Defaults to (10, 10) Controls the size of the matplotlib figure. See [matplotlib documentation](#) for more information.
- **apply_orientation_transform** – Define whether the orientation transform should be applied before plotting.
- **(Optional[str] (save_name) – None)**: Specify a filename for saving the generated figure. By default *None* is provided which will not save a figure.
- **default** – *None*: Specify a filename for saving the generated figure. By default *None* is provided which will not save a figure.

save(*file_location: str, encoding: str = 'utf-8'*)

Can be used to save the shape collection as XML file.

file_location: File path pointing to the XML file.

svg_to_lmd(*file_location, offset=[0, 0], divisor=3, multiplier=60, rotation_matrix=numpy.eye, orientation_transform=None*)

Can be used to save the shape collection as XML file.

Parameters

- **file_location** – File path pointing to the SVG file.
- **orientation_transform** – Will supersede the global transform of the Collection.
- **rotation_matrix** –

2.1.2 Shape

class `lmd.lib.Shape`(*points: numpy.ndarray, well: Optional[str] = None, name: Optional[str] = None, orientation_transform=None*)

Class for creating a single shape object.

from_xml(*root*)

Load a shape from an XML shape node. Used internally for reading LMD generated XML files.

Parameters *root* – XML input node.

to_xml(*id: int, orientation_transform: numpy.ndarray, scale: int*)

Generate XML shape node needed internally for export.

Parameters

- **id** – Sequential identifier of the shape as used in the LMD XML format.

- **orientation_transform** (*np.array*) – Pass orientation_transform which is used if no local orientation transform is set.
- **scale** (*int*) – Scaling factor used to enable higher decimal precision.

Note: If the Shape has a custom orientation_transform defined, the custom orientation_transform is applied at this point. If not, the orientation_transform passed by the parent Collection is used. This highlights an important difference between the Shape and Collection class. The Collection will always has an orientation transform defined and will use *np.eye(2)* by default. The Shape object can have a orientation_transform but can also be set to *None* to use the Collection value.

2.2 lmd.tools

`lmd.tools.get_rotation_matrix(angle: float)`

Returns a rotation matrix for clockwise rotation.

Parameters **angle** (*float*) – Rotation in radian.

Returns Matrix in the shape of (2, 2).

Return type *np.ndarray*

`lmd.tools.glyph(glyph, offset=numpy.array, rotation=0, divisor=10, multiplier=1, **kwargs)`

Get an uncalibrated *lmd.lib.Collection* for a glyph of interest.

Parameters

- **glyph** (*str*) – Single glyph as string.
- **divisor** (*int*) – Parameter which determines the resolution when creating a polygon from a SVG. A larger divisor will lead to fewer datapoints for the glyph. Default value: 10
- **offset** (*np.ndarray*) – Location of the glyph based on the top left corner. Default value: *np.array((0, 0))*
- **multiplier** (*float*) – Scaling parameter for defining the size of the glyph. The default height of a glyph is 10 units. Default value: 1

Returns Uncalibrated Collection which contains the Shapes for the glyph.

Return type *lmd.lib.Collection*

`lmd.tools.glyph_path(glyph)`

Returns the path for a glyph of interest. Raises a *NotImplementedError* if an unknown glyph is requested.

Parameters **glyph** (*str*) – Single glyph as string.

Returns Path for the glyph.

Return type *str*

`lmd.tools.rectangle(width, height, offset=(0, 0), angle=0, rotation_offset=(0, 0))`

Get a *lmd.lib.Shape* for rectangle of choosen dimensions.

Parameters

- **width** (*float*) – Width of the rectangle.
- **offset** (*np.ndarray*) – Location of the rectangle based on the center. Default value: *np.array((0, 0))*

- **angle** (*float*) – Rotation in radian.
- **rotation_offset** (*np.ndarray*) – Location of the center of rotation relative to the center of the rectangle. Default value: `np.array((0, 0))`

Returns Shape which contains the rectangle.

Return type *lmd.lib.Shape*

Example:

```
lmd.tools.text(text, offset=numpy.array, divisor=1, multiplier=1, rotation=0, **kwargs)
```

Get an uncalibrated `lmd.lib.Collection` for a text.

Parameters

- **text** (*str*) – Text as string.
- **divisor** (*int*) – Parameter which determines the resolution when creating a polygon from a SVG. A larger divisor will lead to fewer datapoints for the glyph. Default value: 10
- **offset** (*np.ndarray*) – Location of the text based on the top left corner. Default value: `np.array((0, 0))`
- **multiplier** (*float*) – Scaling parameter for defining the size of the text. The default height of a glyph is 10 units. Default value: 1

Returns Uncalibrated Collection which contains the Shapes for the text.


Return type *lmd.lib.Collection*

Example

```
import numpy as np
from lmd.lib import Collection, Shape
from lmd import tools

calibration = np.array([[0, 0], [0, 100], [100, 50]])
my_first_collection = Collection(calibration_points = calibration)

identifier_1 = tools.text('0456_B2', offset=np.array([30, 40]), rotation = -np.pi/4)
my_first_collection.join(identifier_1)
my_first_collection.plot(calibration = True)
```



PYTHON MODULE INDEX

|

`lmd.tools`, [15](#)

INDEX

A

`add_shape()` (*lmd.lib.Collection method*), 13

C

`calibration_points` (*lmd.lib.Collection attribute*), 13
`Collection` (*class in lmd.lib*), 13

F

`from_xml()` (*lmd.lib.Shape method*), 14

G

`get_rotation_matrix()` (*in module lmd.tools*), 15
`glyph()` (*in module lmd.tools*), 15
`glyph_path()` (*in module lmd.tools*), 15

J

`join()` (*lmd.lib.Collection method*), 13

L

`lmd.tools`
 module, 15
`load()` (*lmd.lib.Collection method*), 13

M

`module`
 lmd.tools, 15

N

`new_shape()` (*lmd.lib.Collection method*), 13

O

`orientation_transform` (*lmd.lib.Collection attribute*),
 13

P

`plot()` (*lmd.lib.Collection method*), 14

R

`rectangle()` (*in module lmd.tools*), 15

S

`save()` (*lmd.lib.Collection method*), 14
`Shape` (*class in lmd.lib*), 14
`shapes` (*lmd.lib.Collection attribute*), 13
`svg_to_lmd()` (*lmd.lib.Collection method*), 14

T

`text()` (*in module lmd.tools*), 16
`to_xml()` (*lmd.lib.Shape method*), 14